Aperta NSE PTA Balloon Dilatation Catheter

EN	English	English	1
DE	German	Deutsch	4
FR	French	Français	8
ES	Spanish	Español	11
IT	Italian	Italiano	14
SV	Swedish	Svenska	17
FI	Finnish	Suomi	20
NL	Dutch	Nederlands	23
PT	Portuguese	Português	26

CS	Czech	Česky	29
SK	Slovak	Slovensky	32
PL	Polish	Polski	35
HU	Hungarian	Magyar	38
RU	Russian	Русский язык	41
RO	Romanian	Limba română	45
EL	Greek	Ελληνικά	48
	Compliance chart		52
	Labeling Symbol Definitions		53

MANUFACTURER:

Goodman Co., Ltd.

5F KDX Nagoya Sakae Building, 4-5-3 Sakae, Naka-ku,
Nagoya, Aichi 460-0008 Japan

E-mail : gm-adm-qa@goodmankk.com

AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY:

Goodman Medical Ireland Ltd

Mervue Business Park Galway H91 H9CK Ireland

EU IMPORTER:

NIPRO MEDICAL EUROPE NV

Blokhuisstraat 42, 2800 Mechelen, Belgium

EU DISTRIBUTOR:

B.Braun Melsungen AG

Carl-Braun-Str.1, Melsungen, 34212 Germany

eIFU:

The latest IFU is available at

<https://www.goodmankk.com/english/ifu/index.php>


18th Feb, 2025/Rev.2/GM142/ L02660002X

ENGLISH

*Instructions for Use (IFU) as per European Regulation (EU) 2017/745

PRODUCT DESCRIPTION

1. Product Diagram

* Products with a working length of 50cm and 90cm do not have depth marks.

<Components>

Catheter clip (Products with a working length of 50cm do not contain the catheter clip.), Rewrap tool / 1 unit

<Material>

Nylon resin, polyethylene resin, polycarbonate, polyethylene oxide, polyvinylpyrrolidone

2. Product Summary

This product is an OTW PTA balloon catheter used for the purpose of dilatation of stenotic lesions during percutaneous transluminal angioplasty.

The balloon expands to the specific OD and length at nominal pressure and has four elements on the outside of the balloon (with elements creating a focused transmission of dilative force, i.e. scoring effect) in order to dilate stenosis that is traditionally considered difficult to expand. This product has a hydrophilic coating applied to the surface.

WARNING

1. Safety and efficacy has not been clearly determined for the use of the product when applied to complex and tortuous lesions in the iliac artery. Careful consideration should be made in regards to the potential for vessel rupture requiring surgical treatment and risk of occlusion of the distal artery from vasculature dissection.
2. Do not exceed balloon expansion beyond the vessel size both proximal and distal to the lesion (damage to vessel may result).
3. Care should be taken when inserting and expanding product at calcified and in-stent restenotic lesions and at location of stiff synthetic grafts (damage to vessel or product may result).
4. Care should be taken to ensure that this product does not become entangled with multiple stents implanted at bifurcated lesions (damage to product may result).
5. Care should be taken when used at the distal location to a drug eluting stent (disruption of newly formed neointimal formation may result).
6. This product should only be manipulated in a completely deflated position (damage to product may result).

CONTRAINDICATIONS

1. Contraindications For Use

This product is sterilized. It is not re-usable (single use only) and must not be re-sterilized. Resterilization and/or reuse could result in infection or degradation of

product characteristics such as balloon size, shaft strength or lubricity and could result in failure of this product during use.

2. Inappropriate Lesions

- (1) Exceeding the stent strut in a bifurcated lesion
- (2) Fractured stent

CLINICAL BENEFITS

The primary purpose of this product is to alleviate the patient's stenosis or ischemic condition of the occluded vessel which causes peripheral complications (claudication and amputation etc.) along with VA trouble (insufficient blood drawing and elevated venous pressure etc.) and results in relief from symptoms. The main purpose of EVT is to maintain long term patency of the lesion following alleviation. This product is a PTA balloon catheter with the same purpose as a conventional balloon or scoring/cutting balloon (generic device group), and is intended to dilate stenotic lesions and occlusions that are determined difficult to dilate by conventional balloons in order to achieve revascularization. The elements of this product create a scoring effect while minimizing dissection with optimal dilatation.

INTENDED USER

Physicians trained in percutaneous intravascular procedures.

INTENDED PURPOSE

This product is intended for dilatation of lesions.

INDICATIONS FOR USE

This product is indicated for percutaneous transluminal angioplasty (PTA) of lesions in peripheral arteries, including iliac, femoral, ilio-femoral, popliteal, and infra-popliteal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This product is not for use in the coronary or neuro-vasculature including carotid arteries.

APPROPRIATE PATIENT SELECTION

Patients with appropriate peripheral arteries and arteriovenous fistula for a balloon diameter in a percutaneous transluminal angioplasty procedure as determined by a physician (including pregnant patients). No clinical data specific for pregnant patients (or the potential for being pregnant) is available.

SPECIFICATIONS

1. Nominal Pressure (NP)

10atm [10×10^2 kPa]

2. Rated Burst Pressure (RBP)

- 1) Balloon diameter 4.0, 5.0, 6.0mm: 20atm [20×10^2 kPa]
- 2) Balloon diameter 7.0, 8.0mm: 14atm [14×10^2 kPa]

* In vitro test results indicate that 99.9% of the balloons do not incur rupture (95% reliability) when expanded up to maximum recommended inflation pressure.

3. The balloon compliance chart is provided at the end of this Instructions for Use

4. Balloon length: 40mm at NP

5. Conformity Of Connector

ISO80369-7

MEDICAL DEVICES REQUIRED FOR THIS PRODUCT

- Guidewire
Maximum guide wire diameter: 0.46mm (0.018inch)
- Sheath introducer

Minimum Introducer Sheath

- 1) Balloon diameter 4.0, 5.0mm: 5F (1.78mm)
- 2) Balloon diameter 6.0, 7.0, 8.0mm: 6F (2.16mm)

- Inflation device with manometer
- Stopcock with luer connection conform to ISO80369-7
- Luer-lock syringe
- Hemostatic valve connector

METHOD OF USE

1. Preparation

- (1) Inspect and confirm each device for correct functionality.
- (2) Remove this product from the packaging (remaining in the protective hoop).
- (3) Carefully remove this product from the protective hoop.
- (4) Carefully remove the balloon protector and the stylet. Confirm that there is no damage to this product.
- (5) Ensure that the entire length of the shaft of the product is saturated in heparinized saline solution.
- (6) Connect the stopcock to this product hub.
- (7) Prime an inflation device with appropriate amount of diluted contrast media (contrast: heparinized saline = 1:1). Connect it to the stopcock and flush the port of the stopcock with contrast media.
- (8) With the tip of the inflation device pointing downwards, remove all residual air within the balloon lumen by continuous application of negative pressure.
- (9) With negative pressure applied to the inflation device, close the stopcock ensuring that negative pressure of the balloon lumen is maintained.

* Do not allow air or liquid to enter the balloon in order to preserve the balloon folding functionality and the elements.

- (10) Attach the syringe to the hub guide wire port, flush the guide wire lumen of the product with heparinized saline, removing all residual air.

2. Insertion Of Product

- (1) Follow the package insert for each medical device and complete preparation of the devices prior to insertion of this product.
- (2) With the balloon completely deflated, insert the guide wire into the distal tip of the product and carefully advance the product to this target lesion.

3. Balloon Inflation

- (1) Confirm position of the balloon in relation to the target lesion and close the hemostatic valve, locking balloon into place.
- (2) While confirming the balloon diameter under fluoroscopy, inflate and deflate the balloon using an inflation device (refer to the compliance chart included in packaging for each balloon size reference diameter at corresponding inflation pressure).
- (3) Perform multiple inflations as deemed necessary.
- (4) Upon completion of dilatation, ensure the balloon is completely deflated, return this product to within the guiding catheter and evaluate improvement to stenotic site via angiography.

4. Removal Of Product

While confirming guide wire positioning, carefully remove the product through the sheath introducer (guiding catheter).

5. Rewrap Tool (reference)

The rewrap tool is used for re-wrapping the balloon folds. When undertaking re-wrapping using the rewrap tool included, follow the procedure described below.

- (1) Insert the stylet through the side of the rewrap tool that does not have a flare.
- (2) With the balloon held in negative pressure, insert the stylet into the guide wire lumen.

- (3) Being careful not to damage product, use fingers to gently roll the balloon wrapping. Looking directly at the fold from the tip, the folds are wrapped in a clock wise direction.
- (4) Carefully insert the balloon body into the rewrap tool. In order to avoid damaging product do not rotate product or tool during insertion.
- (5) While inside the rewrap tool, apply low inflation pressure and slowly reduce pressure.
- (6) Under negative pressure, carefully remove rewrap tool and stylet.
- (7) Visually confirm condition of product (including balloon).

PRECAUTIONS

<Fundamental Considerations>

- (1) This product should be used in those facilities capable of immediate surgical response to any complication with potential negative effects on health or life threatening complication that may arise as a result of its use.
- (2) Use of this product should be undertaken under an appropriate anticoagulation regimen.
- (3) Appropriate balloon sizing and products used in combination should be determined based on the anatomical structure of the patient. Size selection (both diameter and length) should not exceed the proximal and distal sites of the stenotic location, nor exceed lesion length.
- (4) All device manipulation should be carefully performed under angiography while constantly confirming the position of the catheter (damage to vessel and product may result).
- (5) Should any resistance be incurred during usage or any anomaly of the catheter position be observed, immediately cease product manipulation and confirm the cause under angiography (damage to vessel or product may result).
- (6) Should any damage to product be observed, immediately cease use and confirm condition of the elements (damage to vessel may result).
- (7) All product removal from protective packaging and hoop, as well as removal of sheath and stylet, should be performed carefully so as to not place excessive burden on this product (product insertion and functionality may be affected).
- (8) Use a gauze that has been sufficiently dampened with heparinized saline to remove any particles attached to the guide wire when inserting and removing the product (operability of guide wire and product may be affected, damage may result).
- (9) This product should only be used by physicians trained in percutaneous intravascular procedures.
- (10) Do not exceed RBP (damage to product may result).
- (11) All balloon inflations and deflations should be performed under radiology while confirming the entire balloon.
- (12) Do not rotate product (damage may result).
- (13) Do not use products comprising organic solvent, fat-based emulsions or oil (damage to hub and/or elements, as well as deterioration to the hydrophilic coating may result).
- (14) Careful consideration of pharmacological allergies to products used, potential of patient for being pregnant, hemodynamic condition and potential for cardiogenic shock should be undertaken.

1. Precautions Prior To Use

- (1) Use of distal protection is recommended for those patients with potential for distal embolization.
- (2) Refer to attached documents of all pharmaceutical and medical devices used in combination.
- (3) Do not use this product if the packaging or contents are damaged or contaminated, in particular that damage to elements has not incurred.
- (4) All product use should be undertaken in a sterile environment.

2. Precautions During Use

- (1) Do not cause damage to product when using sharp items.
- (2) All product manipulation should be performed with the guide wire advanced from product tip.

- (3) All balloon inflations should be performed with contrast media, ensuring there is no residual air (potential for air embolization may result).
- (4) During balloon inflation confirm condition of inflation device as well as inflation via fluoroscopy. Should any irregularity be observed, immediately cease inflation, deflate balloon and remove the product, confirming potential causes.
- (5) During long periods of use carefully consider thrombus accumulation (product and guide wire functionality may be affected).
- (6) During inflation in a position that does not allow for concentric balloon expansion, care should be taken to ensure the balloon does not move position (damage to vessel may incur).
- (7) Do not overtighten the hemostatic valve, potentially affecting guide wire manipulation, as well as balloon inflation and deflation by impeding flow of contrast media.
- (8) Do not use product if kinking is observed (product may break).
- (9) When using multiple devices in-vivo, ensure that devices do not become entangled. Should any resistance during product delivery be encountered, immediately determine the cause (product damage may result).
- (10) Continuously confirm the condition of the catheter to ensure that there is no damage, that connections are attached and that no leakage of contrast media is determined. Should any abnormality be determined, immediately cease use of this product (complications to the vessel may result).
- (11) Once the product has been removed from the balloon protective sheath, do not re-insert the product (damage to product may result).
- (12) This product has a hydrophilic coating and must remain damp from heparinized saline at all times (operability may be affected and damage may result).

3. Precautions After Use

Dispose of this product as medical waste and take measures to prevent possible spread of infection.

4. Adverse Events

Possible adverse events include, but are not limited to, the following:

- death
- myocardial infarction
- acute closure
- closure, ischemia
- arrhythmia including ventricular fibrillation
- angina
- cerebrovascular disease
- restenosis
- claudication
- necrosis
- quadruple amputation
- embolization (air, tissue, thrombus)
- complications of hemorrhage
- puncture complications
- hematoma
- pain
- vascular spasm
- vessel dissection, perforation, rupture
- arterial dissection
- arteriovenous fistula
- pseudoaneurysm
- infection
- pharmaceutical allergic reaction

STORAGE METHOD, SHELF LIFE AND OTHER

1. Storage Method

- (1) Store the product in a room temperature location not exposed to high temperature and humidity or direct sunlight and take proper precautions to ensure product does not contact water.
- (2) Avoid inclinations, vibrations and impacts (including during transportation) and store in a safe, stable environment.
- (3) Do not store near chemicals or in areas where the product may be exposed to ionizing radiations.

2. Shelf Life

Use this product before the date of "Use By" shown on the package label.

3. Sterilization Method

This product is supplied STERILE using Ethylene Oxide Gas (EtO). Non-pyrogenic.

PACKAGING

1unit/box

PRODUCT WARRANTY DISCLAIMERS AND LIMITATION OF REMEDIES

WITH RESPECT TO ANY AND ALL GOODMAN PRODUCT(S) DEPICTED OR DESCRIBED IN THIS PUBLICATION, THERE ARE NO EXPRESS OR IMPLIED WARRANTIES OF ANY KIND WHATSOEVER, INCLUDING, AND WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GOODMAN CO., LTD AND ITS SUBSIDIARIES ("GOODMAN") SHALL NOT BE LIABLE FOR ANY DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES IN CONNECTION WITH THE USE, RE-USE, OR ANY OTHER ASPECT OF THE RELEVANT PRODUCT(S) EXCEPT AS EXPRESSLY PROVIDED UNDER APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND GOODMAN TO ANY REPRESENTATION OR WARRANTY. ANY PRODUCT OR SERVICE DESCRIPTION OR SPECIFICATIONS CONTAINED IN ANY GOODMAN PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE RELEVANT PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS OR OTHER WARRANTIES.

Any serious incidents shall be reported to the manufacturer and the competent authority of your State.

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