Goodnet[®] Version7

DICOM Conformance Statement

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2024. 3

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表 1-2 Supported Transfer Syntaxes

Transfer Syntax Set	Transfer Syntax Name	Transfer Syntax UID	DICOM Web Service Bulkdata Media Type
Lossless Compressed Transfer Syntax Set (LL)	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1])	1.2.840.10008.1.2.4.70	Image/jpeg
	RLE Lossless	1.2.840.10008.1.2.5	Image/x-dicom-rle
Lossy Compressed Transfer Syntax Set (L)	JPEG Baseline (Process 1)	1.2.840.10008.1.2.4.50	Image/jpeg
	JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	Image/jpeg
Non-Image Transfer Syntax Set (NI)	Implicit VR Little Endian	1.2.840.10008.1.2	
	Explicit VR Little Endian	1.2.840.10008.1.2.1	
	Explicit VR Big Endian	1.2.840.10008.1.2.2	
Uncompressed Transfer Syntax Set (U)	Implicit VR Little Endian	1.2.840.10008.1.2	
	Explicit VR Little Endian	1.2.840.10008.1.2.1	
	Explicit VR Big Endian (Retired)	1.2.840.10008.1.2.2	

1.1.1. Structured Reporting Root Template IDs

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表 1-5 Workflow Management SOP Classes

SOP Classes		Transfer Syntax		SCU	SCP
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian	1.2.840.10008.1.2	※Y	N
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

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Service	Transaction	Resource	User Agent	Origin Server
		Bulkdata	N	N
	Search (QIDO-RS)	All Studies	※Y	N
		Study	※Y	N
		Study's Series	※Y	N
		Study's Instances	N	N
		All Series	N	N
		Series	N	N
		Series Instances	※Y	N
		All Instances	N	N
		Instance	N	N
		Store (STOW-RS)	All Study	N
	Study		N	N
	Bulkdata		N	N

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N/A

1.7. Specific Character Sets

表 1-11 Supported Specific Character Sets

Defined Term	IANA	Description
Single-Byte Character Sets without Code Extensions		
ISO_IR 6	ISO-646 or US-ASCII	Default Repertoire
ISO_IR 13	JIS_X0201	Japanese JIS X 0201
Single-Byte Character Sets with Code Extensions		
ISO 2022 IR 6		Default Repertoire
ISO 2022 IR 13		Japanese JIS X 0201
Multi-Byte Character Sets with Code Extensions		
ISO 2022 IR 87	ISO-2022-JP	Japanese JIS X 0208
ISO 2022 IR 159	ISO-2022-JP-2	Japanese JIS_X0212-1990
Multi-Byte Character Sets without Code Extensions		
ISO_IR 192	UTF-8	Unicode in UTF-8

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1.	概要	2
1.1.	コンテンツと転送.....	3
1.1.1.	Structured Reporting Root Template IDs.....	6
1.2.	DIMSE Services	7
1.2.1.	Verification.....	7
1.2.2.	Storage	7
1.2.3.	Workflow Management.....	7
1.2.4.	Query/Retrieve	8
1.2.5.	Printing - N/A	8
1.3.	DIMSE Web Services	8
1.3.1.	URI Service (WADO-URI).....	8
1.3.2.	Studies Service.....	9
1.3.3.	Worklist Service - N/A.....	10
1.3.4.	Non-Patient Service - N/A.....	10
1.3.5.	Storage Commitment Service - N/A.....	10
1.4.	Media Services	11
1.5.	Real Time Video Services - N/A	11
1.6.	De-identification Profile - N/A.....	11
1.7.	Specific Character Sets.....	12
2.	目次	13
3.	まえがき	17
3.1.	改訂履歴	17
3.2.	対象者.....	17
3.3.	備考	18
3.4.	用語と定義.....	19
3.5.	略語	21
3.6.	参考文献	23
4.	実装モデル.....	24
4.1.	アプリケーションデータフロー図.....	24
4.1.1.	A Eの機能定義	25

5.	サービスと相互運用性の説明.....	28
5.1.	サービスとアプリケーション・エンティティのマッピング.....	28
5.2.	DIMSEサービス.....	29
5.2.1.	Basic Worklist Management Service.....	29
5.2.2.	Modality Performed Procedure Step Service - N/A.....	31
5.2.3.	Unified Worklist and Procedure Step Service - N/A.....	31
5.2.4.	Instance Availability Notification SOP Class - N/A.....	31
5.2.5.	Storage Service.....	31
5.2.6.	Storage Commitment Service.....	34
5.2.7.	Query/Retrieve Service Class.....	35
5.2.8.	Print Management Service - N/A.....	41
5.3.	Supported DICOM Web Services.....	42
5.3.1.	URI Web Service (WADO-URI).....	42
5.3.2.	Studies Web Service.....	43
5.3.3.	Worklist Web Service - N/A.....	46
5.3.4.	Non-Patient Instance Web Service - N/A.....	47
5.3.5.	Notification Web Service - N/A.....	47
5.4.	Media Services.....	47
5.4.1.	File Set Creator (FSC).....	47
5.4.2.	File Set Reader (FSR).....	47
5.4.3.	File Set Updater (FSU) - N/A.....	47
5.5.	Real Time Video Services - N/A.....	47
5.6.	Cross Service Considerations - N/A.....	47
5.7.	Specific Character Set.....	47
6.	Configuration.....	49
6.1.	General Configuration Parameters.....	49
6.2.	Configuration of DIMSE Services.....	50
6.2.1.	Basic Worklist Management Service Configuration.....	50
6.2.2.	Modality Performed Procedure Step Service Configuration - N/A -.....	51
6.2.3.	Unified Worklist and Procedure Step Service Configuration - N/A -.....	51
6.2.4.	Instance Availability Notification Service Configuration - N/A -.....	51
6.2.5.	Storage Service Configuration.....	51
6.2.6.	Storage Commitment Service Configuration.....	52
6.2.7.	Query/Retrieve Service Configuration.....	53
6.2.8.	Print Management Service Configuration - N/A -.....	55

6.3.	Configuration of DICOM Web Services.....	55
6.3.1.	URI Web Service Configuration.....	55
6.3.2.	Studies Web Service Configuration.....	56
6.3.3.	Worklist Web Service Configuration - N/A	56
6.3.4.	Non-Patient Instances (NPI) Web Service Configuration - N/A.....	56
6.3.5.	Storage Commitment Service Configuration - N/A.....	56
6.4.	Configuration of Media Storage Service	57
6.5.	Configuration of Real Time Video Service - N/A	57
6.6.	Configuration of Audit Trail – Syslog - N/A.....	57
7.	ネットワークおよびメディア交換の詳細.....	58
7.1.	General.....	58
7.1.1.	General Association Parameters - N/A.....	59
7.2.	仕様	59
7.2.1.	DICOM受信サービス Application Entity.....	59
7.2.2.	DICOM自動転送サービス Application Entity	63
7.2.3.	DICOM Q/R SCP サービス Application Entity	65
7.2.4.	DICOM Q/R SCUサービス Application Entity.....	68
7.2.5.	DICOM画像送信サービス Application Entity	71
7.2.6.	DICOM MWL SCUサービス Application Entity.....	73
7.2.7.	DICOM Webクライアント Application Entity.....	75
7.2.8.	メディア保存サービス Application Entity	75
7.2.9.	メディア登録サービス Application Entity	76
7.3.	Status Code.....	77
7.3.1.	一般的なAEコミュニケーションと失敗の挙動と対応.....	77
7.3.2.	DIMSE Services.....	78
7.3.3.	DICOM Web Service.....	87
8.	セキュリティ	91
8.1.	はじめに	91
8.2.	外部ネットワーク要件.....	91
8.3.	TCP Port Configuration.....	91
8.4.	DICOM Security Profiles Support – N/A	91
8.5.	User Identity Negotiation Support – N/A	91
8.6.	Web Services Security Features – N/A.....	91
8.7.	Other Security Features – N/A.....	91

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AE	Application Entity
AET	Application Entity Title
CAD	Computer Aided Detection
CDA	Clinical Document Architecture
CID	Context Identifier
DCS	DICOM Conformance Statement
DHCP	Dynamic Host Configuration Protocol
DICOM	Digital Imaging and Communications in Medicine
ELE	Explicit VR Little Endian
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
GUI	Graphical User Interface
HTTP	Hyper Text Transfer Protocol
IANA	Internet Assigned Numbers Authority
IHE	Integrating the Healthcare Enterprise
ILE	Implicit VR Little Endian
IOD	Information Object Definition
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standardization
JSON	JavaScript Object Notation

MPPS	Modality Performed Procedure Step
MWL	Modality Worklist
NEMA	National Electrical Manufacturers Association
NTP	Network Time Protocol
OID	Object Identifier
OS	Origin Server
PACS	医用画像保存通信システム (Picture Archiving and Communication System)
PDU	Protocol Data Unit
PHI	Protected Health Information
PPS	Performed Procedure Step
QIDO-RS	Query based on ID for DICOM Objects by RESTful Services
REST	Representational State Transfer
RIS	放射線科情報システム (Radiology Information System)
RTV	Real Time Video
SCP	Service Class Provider
SCU	Service Class User
SDP	Service Description Protocol
SOP	Service-Object Pair
SPS	Scheduled Procedure Step
SR	Structured Reporting
STOW-RS	STore Over the Web by RESTful Services
TCP/IP	Transmission Control Protocol/Internet Protocol
TID	Template Identifier
UA	User Agent
UI	User Interface
UID	Unique Identifier
UL	Upper Layer
UPS	Unified Procedure Step
UPS-RS	Unified Procedure Step by RESTful Services
URL	Uniform Resource Locator
VR	Value Representation
WADO-RS	Web Access to DICOM Objects by RESTful Services
WADO-URI	Web Access to DICOM Objects by URI

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- [1] National Electrical Manufacturers Association (NEMA), Rosslyn, VA USA. PS3 / ISO 12052 Digital Imaging and Communications in Medicine (DICOM) Standard. <http://www.dicomstandard.org>.
- [2] Integrating the Healthcare Enterprise (IHE). IHE Technical Framework. http://www.ihe.net/Resources/technical_frameworks/

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表 5-2 Supported C-FIND Query Parameters for Modality Worklist - SCU

Attribute Name	Tag	Matching Type	Query Value Source	Value	Display on UI	Comments
Schedule Procedure Step						
Scheduled Procedure Step Sequence	(0040,0100)	SEQUENCE			N	
>Scheduled Station AE Title	(0040,0001)	SINGLE_VALUE	CONFIGURATION		D	ワークリストを取得するシステムの AE タイトル
>Scheduled Procedure Step Start Date	(0040,0002)	RANGE	CONFIGURATION		D	
>Scheduled Procedure Step Start Time	(0040,0003)	UNIVERSAL			D	
>Scheduled Procedure Physician's Name	(0040,0006)	UNIVERSAL			D	
>Scheduled Procedure Step Description	(0040,0007)	UNIVERSAL			D	
>Scheduled Protocol Code Sequence	(0040,0008)	UNIVERSAL			N	
>>Code Value	(0008,0100)	UNIVERSAL			D	
>>Coding Scheme Designer	(0008,0103)	UNIVERSAL			D	
>>Coding Scheme Version	(0008,0104)	UNIVERSAL			D	
>>Code Meaning	(0008,0104)	UNIVERSAL			D	
>Scheduled Procedure Step ID	(0040,0009)	UNIVERSAL			D	
>Modality	(0008,0060)	SINGLE VALUE			D	
Requested Procedure ID	(0040,1001)	UNIVERSAL			D	
Requested Procedure Code Sequence	(0032,1064)	UNIVERSAL			N	
Study Instance UID	(0040,000D)	UNIVERSAL			D	
Requested Procedure	(0032,1060)	UNIVERSAL			D	
Referring Physician	(0032,1032)	UNIVERSAL			D	
Requesting Service	(0032,1033)	UNIVERSAL			D	
Accession Number	(0008,0050)	UNIVERSAL			D	
Patient's Name	(0010,0010)	WILDDCARD			D	
Patient ID	(0010,0020)	UNIVERSAL			D	
Patient's Birth Date	(0010,0030)	UNIVERSAL			D	
Patient's Sex	(0010,0040)	UNIVERSAL			D	
Patient's Weight	(0010,1030)	UNIVERSAL			D	

5.2.1.2. SCP of Modality Worklist Information Model - FIND SOP Class - N/A

N/A

5.2.2.Modality Performed Procedure Step Service - N/A

N/A

5.2.3.Unified Worklist and Procedure Step Service - N/A

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表 5-3 Transcoding of Transfer Syntaxes

Stored Transfer Syntax	Sent Transfer Syntax						
	Implicit VR Little Endian	Explicit VR Little Endian	Explicit VR Big Endian	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14)	JPEG Baseline (Process 1)	JPEG Extended (Process 2 & 4)	RLE Lossless
Implicit VR Little Endian		SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED
Explicit VR Little Endian	NOT_SUPPORTED		NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED
Explicit VR Big Endian	NOT_SUPPORTED	SUPPORTED		NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED
JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14)	NOT_SUPPORTED	SUPPORTED	NOT_SUPPORTED		NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED
JPEG Baseline (Process 1)	NOT_SUPPORTED	SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED		NOT_SUPPORTED	NOT_SUPPORTED
JPEG Extended (Process 2 & 4)	NOT_SUPPORTED	SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED		NOT_SUPPORTED
RLE Lossless	NOT_SUPPORTED	SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	NOT_SUPPORTED	

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表 5-5 Levels of Conformance

Levels of Conformance	2
Level of Digital Signature Support	1

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表 5-6 Attribute Coercion by Storage SCP

Attribute Name	Tag	SOP Class UID	Type of Change	New Value	Condition	Comments
Patient's Name	(0010,0020)	ALL	MODIFIED	Order Value	EXTERNAL	
Patient's BirthDate	(0010,0030)	ALL	MODIFIED	Order Value	EXTERNAL	
Patient's Sex	(0010,0040)	ALL	MODIFIED	Order Value	EXTERNAL	

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表 5-7 Display and Processing Limitations for Storage SCP

Limitation Case			Effect	Comments
Attribute Name	Tag	Value		
All SOP Class				
Bits Allocated	(0028,0100)	16	LD	Photometric Interpretation(0028,0004)=RGBは未対応

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5.2.6.Storage Commitment Service

5.2.6.1. SCU of the Storage Commitment SOP Class - N/A

N/A

5.2.6.2. SCP of the Storage Commitment SOP Class

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表 5-10 Failure Conditions on Storage Commitment SCP

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5.2.7. Query/Retrieve Service Class

5.2.7.1. SCU of the StudyRoot Q/R Information Model - FIND SOP Class

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表 5-11 Supported C-FIND Attribute Matching for Study Root Q/R Model - SCU

Attribute Name	Tag	Matching Type	Query Value Source	Value	Display on UI	Comments
Study Level						
Specific Character Set	(0008,0005)	SINGLE_VALUE	GENERATED		N	使用する文字によっては、この属性が入ることがあります。
Study Date	(0008,0020)	RANGE	USER		D	
Study Time	(0008,0030)	UNIVERSAL	EMPTY		D	
Accession Number	(0008,0050)	WILDCARD	USER		D	
Query Retrieve Level	(0008,0052)	SINGLE_VALUE	GENERATED	STUDY	N	
Modalities in Study	(0008,0061)	UNIVERSAL	EMPTY		D	
Patient's Name	(0010,0010)	WILDCARD	USER		D	
Patient ID	(0010,0020)	WILDCARD	USER		D	
Patient Birth Date	(0010,0030)	RANGE	USER		D	
Patient Sex	(0010,0040)	UNIVERSAL	EMPTY		D	
Study Instance UID	(0020,000D)	UNIVERSAL	EMPTY		D	
Study ID	(0020,0010)	UNIVERSAL	EMPTY		D	
Series Level						
Series Date	(0008,0021)	UNIVERSAL	EMPTY		C	
Series Time	(0008,0031)	UNIVERSAL	EMPTY		C	
Query Retrieve Level	(0008,0052)	SINGLE_VALUE	GENERATED	SERIES	N	
Modality	(0008,0060)	SINGLE_VALUE	USER		D	Option
Study Instance UID	(0020,000D)	SINGLE_VALUE	FIXED		D	Unique Key
Series Instance UID	(0020,000E)	UNIVERSAL	EMPTY		D	
Series Number	(0020,0011)	UNIVERSAL	EMPTY		D	
Instance Level						
SOP Instance UID	(0008,0018)	UNIVERSAL	EMPTY		D	
Query Retrieve Level	(0008,0052)	SINGLE_VALUE	GENERATED	IMAGE	N	

Attribute Name	Tag	Matching Type	Query Value Source	Value	Display on UI	Comments
Study Instance UID	(0020,000D)	SINGLE_VALUE	FIXED		D	Unique Key
Series Instance UID	(0020,000E)	SINGLE_VALUE	FIXED		D	Unique Key
Instance Number	(0020,0013)	UNIVERSAL	EMPTY		D	
Number of Frames	(0028,0008)	UNIVERSAL	EMPTY		C	

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5.2.7.2. SCU of the PatientRoot Q/R Information Model - FIND SOP Class - N/A

N/A

5.2.7.3. SCU of the StudyRoot Q/R Information Model - MOVE SOP Class

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表 5-12 Supported C-FIND Attribute Matching for Study Root Q/R Model - SCP

Attribute Name	Tag	Matching Type	Comments
Study Level			
Specific Character Set	(0008,0005)	SINGLE_VALUE	
Study Date	(0008,0020)	RANGE	
Study Time	(0008,0030)	RANGE	
Accession Number	(0008,0050)	WILDCARD	
Modalities in Study	(0008,0061)	SINGLE_VALUE	
SOP Classes in Study	(0008,0062)	UNIVERSAL	
Referring Physician's Name	(0008,0090)	WILDCARD	
Study Description	(0008,1030)	WILDCARD	
Patient's Name	(0010,0010)	WILDCARD	
Patient ID	(0010,0020)	SINGLE_VALUE	
Patient's Birth Date	(0010,0030)	RANGE	
Patient's Birth Time	(0010,0032)	RANGE	
Patient's Sex	(0010,0040)	SINGLE_VALUE	
Patient's Age	(0010,1010)	SINGLE_VALUE	
Patient's Size	(0010,1020)	SINGLE_VALUE	
Patient's Weight	(0010,1030)	SINGLE_VALUE	
Patient's Comments	(0010,4000)	WILDCARD	
Study Instance UID	(0020,000D)	UID	
Study ID	(0020,0010)	WILDCARD	
Number of Study Related Series	(0020,1206)	UNIVERSAL	
Number of Study Related Instances	(0020,1208)	UNIVERSAL	
Series Level			
Specific Character Set	(0008,0005)	SINGLE_VALUE	
Modality	(0008,0060)	SINGLE_VALUE	
Series Description	(0008,103E)	WILDCARD	
Study Instance UID	(0020,000D)	SINGLE_VALUE	Unique Key
Series Instance UID	(0020,000E)	UID	

Attribute Name	Tag	Matching Type	Comments
Number of Series Related Instances	(0020,1209)	UNIVERSAL	
Instance Level			
SOP Class UID	(0008,0016)	UID	
SOP Instance UID	(0008,0018)	UID	
Study Instance UID	(0020,000D)	SINGLE_VALUE	Unique Key
Series Instance UID	(0020,000E)	SINGLE_VALUE	Unique Key
Instance Number	(0020,0013)	SINGLE_VALUE	
Number of Frames	(0028,0008)	SINGLE_VALUE	
Rows	(0028,0010)	SINGLE_VALUE	
Columns	(0028,0011)	SINGLE_VALUE	
Bits Allocated	(0028,0100)	SINGLE_VALUE	

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表 5-13 Supported C-FIND Attribute Matching for PatientRoot Q/R Model - SCP

Attribute Name	Tag	Matching Type	Comments
Patient Level			
Specific Character Set	(0008,0005)	SINGLE_VALUE	
Patient's Name	(0010,0010)	WILDCARD	
Patient ID	(0010,0020)	WILDCARD	
Patient's Birth Date	(0010,0030)	RANGE	
Patient's Birth Time	(0010,0032)	RANGE	
Patient's Sex	(0010,0040)	SINGLE_VALUE	
Patient's Comments	(0010,4000)	WILDCARD	
Number of Patient Related Studies	(0020,1200)	UNIVERSAL	
Number of Patient Related	(0020,1202)	UNIVERSAL	

Attribute Name	Tag	Matching Type	Comments
Series			
Number of Patient Related Instances	(0020,1204)	UNIVERSAL	
Study Level			
Specific Character Set	(0008,0005)	SINGLE_VALUE	
Study Date	(0008,0020)	RANGE	
Study Time	(0008,0030)	RANGE	
Accession Number	(0008,0050)	WILDCARD	
Modalities in Study	(0008,0061)	SINGLE_VALUE	
SOP Classes in Study	(0008,0062)	UNIVERSAL	
Referring Physician's Name	(0008,0090)	WILDCARD	
Study Description	(0008,1030)	WILDCARD	
Patient ID	(0010,0020)	SINGLE_VALUE	Unique Key
Patient's Age	(0010,1010)	SINGLE_VALUE	
Patient's Size	(0010,1020)	SINGLE_VALUE	
Patient's Weight	(0010,1030)	SINGLE_VALUE	
Study Instance UID	(0020,000D)	UID	
Study ID	(0020,0010)	WILDCARD	
Number of Study Related Series	(0020,1206)	UNIVERSAL	
Number of Study Related Instances	(0020,1208)	UNIVERSAL	
Series Level			
Specific Character Set	(0008,0005)	SINGLE_VALUE	
Modality	(0008,0060)	SINGLE_VALUE	
Series Description	(0008,103E)	WILDCARD	
Patient ID	(0010,0020)	SINGLE_VALUE	Unique Key
Study Instance UID	(0020,000D)	SINGLE_VALUE	Unique Key
Series Instance UID	(0020,000E)	UID	
Number of Series Related Instances	(0020,1209)	UNIVERSAL	
Instance Level			
SOP Class UID	(0008,0016)	UID	
SOP Instance UID	(0008,0018)	UID	
Patient ID	(0010,0020)	SINGLE_VALUE	Unique Key
Study Instance UID	(0020,000D)	SINGLE_VALUE	Unique Key
Series Instance UID	(0020,000E)	SINGLE_VALUE	Unique Key
Instance Number	(0020,0013)	SINGLE_VALUE	
Number of Frames	(0028,0008)	SINGLE_VALUE	
Rows	(0028,0010)	SINGLE_VALUE	
Columns	(0028,0011)	SINGLE_VALUE	
Bits Allocated	(0028,0100)	SINGLE_VALUE	

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Header Field	Supported Values	Comments
Accept	image/jpeg	

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表 5-17 Resources Retrieve Transaction – User Agent

Resource	Comments
DICOM Instance Resources	
Study Instances	
Series Instances	
Individual Instance	
DICOM Metadata Resources	
Instance Metadata	

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表 5-19 Header Fields for Retrieve Transaction - User Agent

Header Field	Supported Values	Comments
Instance resource		
Accept	multipart/related; type="application/dicom"	
Metadata resource		
Accept	application/dicom+json	Default
	multipart/related; type="application/dicom+xml"	Option
Cache-control	no-cache	

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表 5-20 Resources Search Transaction - User Agent

Resource	Comments
	See resource path in Table 10.6.1-1 in PS3.18
All Studies	
Study's Series	
Study Series's Instances	

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表 5-23 Header Fields for Search Transaction - User Agent

Header Field	Supported Values	Comments
Accept	application/dicom+json	Default
	multipart/related; type="application/dicom+xml"	Option
Cache-control	no-cache	

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6.2.1. Basic Worklist Management Service Configuration

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6.2.2.Modality Performed Procedure Step Service Configuration - N/A -

N/A

6.2.3.Unified Worklist and Procedure Step Service Configuration - N/A -

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6.2.6.Storage Commitment Service Configuration

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6.2.8. Print Management Service Configuration - N/A -

N/A

6.3. Configuration of DICOM Web Services

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6.3.2.2. Store Transaction (STOW-RS) Configuration - N/A

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6.3.4. Non-Patient Instances (NPI) Web Service Configuration - N/A

N/A

6.3.5. Storage Commitment Service Configuration - N/A

N/A

6.4. Configuration of Media Storage Service

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6.5. Configuration of Real Time Video Service - N/A

N/A

6.6. Configuration of Audit Trail – Syslog - N/A

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表 7-10 Transfer Syntax Selection Preference Order - Image SOP Classes for SCP of the Storage

Preference Order	Transfer Syntax	UID	Comments
1	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection	1.2.840.10008.1.2.4.70	
2	JPEG Baseline (Process 1)	1.2.840.10008.1.2.4.50	
3	JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	
4	Explicit Little-Endian	1.2.840.10008.1.2.1	
5	Explicit Big-Endian	1.2.840.10008.1.2.2	
6	Implicit Little-Endian	1.2.840.10008.1.2	
7	RLE Lossless	1.2.840.10008.1.2.5	

表 7-11 Transfer Syntax Selection Preference Order - Non-Image SOP Classes for SCP of the Storage

Preference Order	Transfer Syntax	UID	Comments
1	Explicit VR Little-Endian	1.2.840.10008.1.2.1	
2	Explicit Big-Endian	1.2.840.10008.1.2.2	
3	Implicit Little-Endian	1.2.840.10008.1.2	

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	Name	Value
Networking Service	Application Context Name	1.2.840.10008.3.1.1.1
	Implementation Class UID	2.16.840.1
	Implementation Version Name	MergeCOM3_460
	Maximum PDU Length	Default: 64234
	ARTIM Timeout	Default: 30s
	Maximum number of simultaneous Associations as Association Initiator	1

	Maximum number of simultaneous Associations as Association Acceptor	1
	Maximum number of outstanding asynchronous Transactions	非同期の操作はサポートしません。

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表 7-14 Transfer Syntax Selection Preference Order – Modality Worklist Management SCU

Preference Order	Transfer Syntax	UID	Comments
1	Explicit VR Little-Endian Transfer Syntax	1.2.840.10008.1.2.1	
2	Implicit VR Little-Endian Transfer Syntax	1.2.840.10008.1.2	
3	Explicit VR Big-Endian Transfer Syntax	1.2.840.10008.1.2.2	

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表 7-17 DICOM Communication Failure Behavior as Association Initiator

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表 7-18 DICOM Communication Failure Handling as Association Acceptor

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7.3.2. DIMSE Services

7.3.2.1. Basic Worklist Management Service

7.3.2.1.1. SCU of the Modality Worklist Information Model Find SOP Class - C-FIND

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表 7-19 Status Codes for C-FIND of the Modality Worklist Information Model SOP Class - SCU

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7.3.2.1.2. SCP of the Modality Worklist Information Model Find SOP Class - C-FIND – N/A

N/A

7.3.2.2. Modality Performed Procedure Step Service – N/A

N/A

7.3.2.3. Unified Worklist and Procedure Step Service – N/A

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7.3.2.6. Storage Commitment Service

7.3.2.6.1. SCU of the Storage Commitment Push Model SOP Class - N-ACTION - N/A

N/A

7.3.2.6.2. SCU of the Storage Commitment Push Model SOP Class - N-EVENT-REPORT - N/A

N/A

7.3.2.6.3. SCP of the Storage Commitment Push Model SOP Class - N-ACTION

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表 7-23 Status Codes for N-ACTION of the Storage Commitment Push Model SOP Class - SCP

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7.3.2.6.4. SCP of the Storage Commitment Push Model SOP Class - N-EVENT-REPORT

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表 7-24 Status Codes for N-EVENT-REPORT of the Storage Commitment Push Model SOP Class - SCP

Service Status	Further Meaning	Status Code	Behavior
----------------	-----------------	-------------	----------

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7.3.3.1.1. Common Transaction as Origin Server – N/A

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表 7-31 Status Codes of User Agent for All Transactions

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表 7-32 Status Codes of User Agent for URI Service

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7.3.3.3. Studies Web Service

7.3.3.3.1. Retrieve Transaction as Origin Server - N/A

N/A

7.3.3.3.2. Retrieve Transaction as User Agent

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7.3.3.3.3. Store Transaction as Origin Server - N/A

N/A

7.3.3.3.4. Store Transaction as User Agent - N/A

N/A

7.3.3.3.5. Search Transaction as Origin Server - N/A

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8.4. DICOM Security Profiles Support – N/A

N/A

8.5. User Identity Negotiation Support – N/A

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Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Image Type	(0008,0008)	FIXED/ SRC_INSTANCE	ALWAYS	CONDITIO NAL	患者検査特性を「SECONDARY」それ以外の項目は、元の値		
SOP Class UID	(0008,0016)	FIXED	ALWAYS	ALWAYS	1.2.840.10008.5.1.4.1.1.7		Secondaly Capture Image Storage
SOP Instance UID	(0008,0018)	GENERATED	ALWAYS	ALWAYS	1.2.392.200082.1.4...		
Conversion Type	(0008,0064)	FIXED	ALWAYS	ALWAYS	WSD		
Date of Secondary Capture	(0018,1012)	GENERATED	ALWAYS	ALWAYS	属性を変更した日付		
Time of Secondary Capture	(0018,1014)	GENERATED	ALWAYS	ALWAYS	属性を変更した時刻		
Secondary Capture Device Manufacturer	(0018,1016)	FIXED	ALWAYS	ALWAYS	GOODMAN		
Secondary Capture Device Software Versions	(0018,1019)	FIXED	ALWAYS	ALWAYS	GN_SCU_7.XX		XX は、ソフトウェアの管理番号
Instance Number	(0020,0013)	FIXED	ALWAYS	ALWAYS	1		
Number of Frames	(0028,0008)	FIXED	ALWAYS	ALWAYS	1		

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Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Media Storage SOP Instance UID	(0002,0003)	GENERATED	ALWAYS	ALWAYS	1.2.392.200082.33 3.302. ...		
Implementation Class UID	(0002,0012)	FIXED	ALWAYS	ALWAYS	1.2.392.200082.1		
Implementation Version Name	(0002,0013)	FIXED	ALWAYS	ALWAYS	GN_PMC7.XX		XX は、ソフトウェアの管理番号
<i>File Set Identification Module</i>							
File Set ID	(0004,1130)	FIXED	ALWAYS	ALWAYS	GOODMAN		
Offset of the First Directory Record of the Root Directory Entity	(0004,1200)	GENERATED	ALWAYS	ALWAYS			
Offset of the Last Directory Record of the Root Directory Entity	(0004,1202)	GENERATED	ALWAYS	ALWAYS			
File-set Consistency Flag	(0004,1212)	GENERATED	ALWAYS	ALWAYS			
Directory Record Sequence	(0004,1220)	GENERATED	ALWAYS	ALWAYS			
Offset of the Next Directory Record	(0004,1400)	GENERATED	ALWAYS	ALWAYS			
Record In-use Flag	(0004,1410)	GENERATED	ALWAYS	ALWAYS			
Offset of Referenced Lower-Level Directory Entity	(0004,1420)	GENERATED	ALWAYS	ALWAYS			
Directory Record Type	(0004,1430)	GENERATED	ALWAYS	ALWAYS			
Referenced File ID	(0004,1500)	GENERATED	ALWAYS	ALWAYS			
Referenced SOP Class UID in File	(0004,1511)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Referenced Transfer Syntax UID in File	(0004,1512)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
<i>Patient Keys</i>							
Specific Character Set	(0008,0005)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Patient's Name	(0010,0010)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Patient ID	(0010,0020)	SRC_INSTANCE / FIXED	SRC_COPY	SRC_COPY	空の場合は「UNKOWN」をセット		
<i>Study Keys</i>							
Study Date	(0008,0020)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study Time	(0008,0030)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study Description	(0008,1030)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study Instance UID	(0020,000d)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study ID	(0020,0010)	SRC_INSTANCE / FIXED	SRC_COPY	SRC_COPY	空の場合は「UNKOWN」をセット		
Accession Number	(0008,0050)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
<i>Series Keys</i>							
Specific Character Set	(0008,0005)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Modality	(0008,0060)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Series Instance UID	(0020,000E)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Series Number	(0020,0011)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
<i>Image Keys</i>							
Specific Character Set	(0008,0005)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Instance Number	(0020,0013)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Samples per Pixel	(0028,0002)	FIXED	ALWAYS	ALWAYS	1		

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Photometric Interpretation	(0028,0004)	FIXED	ALWAYS	CONDITIONAL	「MONOCHROME2」 「PALETTE COLOR」 のいずれかをセット	DICOM を元に判定	
Rows	(0028,0010)	FIXED	ALWAYS	ALWAYS	128		
Columns	(0028,0011)	FIXED	ALWAYS	ALWAYS	128		
Bits Allocated	(0028,0100)	FIXED	ALWAYS	ALWAYS	8		
Bits Stored	(0028,0101)	FIXED	ALWAYS	ALWAYS	8		
High Bit	(0028,0102)	FIXED	ALWAYS	ALWAYS	7		
Pixel Representation	(0028,0103)	FIXED	ALWAYS	ALWAYS	0		
<i>SR Document Keys</i>							
Specific Character Set	(0008,0005)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Instance Number	(0020,0013)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Completion Flag	(0040, A491)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Verification Flag	(0040, A493)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Content Date	(0008,0023)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Content Time	(0008,0033)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Verification Date Time	(0040,A030)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Concept Name Code Sequence	(0040,A043)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Code Value	(0008,1000)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Coding Scheme Designator	(0008,1002)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Coding Scheme Version	(0008,1003)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Code Meaning	(0008,1004)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

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表 A-3-2 General Study Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Study Date	(0008,0020)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Study Time	(0008,0030)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Accession Number	(0008,0050)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study Description	(0008,1030)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study Instance UID	(0020,000d)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Study ID	(0020,0010)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

表 A-3-3 Patient Study Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Patient's Age	(0010,1010)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

表 A-3-4 General Series Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Series Date	(0008,0021)	SRC_INSTANCE /GENERATED	SRC_COPY /ALWAYS	CONDITIO NAL			通常は元の値。 SeriesInstanceUIDを 生成した場合は DICOM 作成日付
Series Time	(0008,0031)	SRC_INSTANCE /GENERATED	SRC_COPY /ALWAYS	CONDITIO NAL			通常は元の値。 SeriesInstanceUIDを 生成した場合は DICOM 作成時刻
Modality	(0008,0060)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Series Instance UID	(0020,000e)	SRC_INSTANCE /GENERATED	ALWAYS	CONDITIO NAL	1.2.392.200082.A. 304.B		通常は元の値。 SeriesInstanceUIDを 生成した場合は、 VALUEの値。
Series Number	(0020,0011)	SRC_INSTANCE /GENERATED	SRC_COPY /ALWAYS	CONDITIO NAL			通常は元の値。 SeriesInstanceUIDを 生成した場合は、元の の値+10000。

表 A-3-5 General Equipment Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Manufacturer	(0008,0070)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Manufacturer Model Name	(0008,1090)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

表 A-3-6 SC Equipment Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Conversion Type	(0008,0064)	FIXED	ALWAYS	ALWAYS	SYN		Synthetic Image
Secondary Capture Device Manufacturers	(0018,1016)	FIXED	ALWAYS	ALWAYS	GOODMAN Co., LTD.		
Secondary Capture Device Manufacturers Model Name	(0018,1018)	FIXED	ALWAYS	ALWAYS	G-Joint		下肢結合を有するアプリケーション名
Secondary Capture Device Software Versions	(0018,1019)	GENERATED	ALWAYS	ALWAYS	X.X.X.X		XはG-Jointのバージョン番号

表 A-3-7 General Image Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Image Type	(0008,0008)	FIXED	ALWAYS	ALWAYS	DERIVED¥SECONDARY		
Acquisition Date	(0008,0022)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Content Date	(0008,0023)	GENERATED	ALWAYS	ALWAYS	DICOM を生成した日付		
Acquisition Time	(0008,0032)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Content Time	(0008,0033)	GENERATED	ALWAYS	ALWAYS	DICOM を生成した時刻		
Instance Number	(0020,0013)	GENERATED	ALWAYS	ALWAYS	元の値に 10000 を足した値		

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Image Comments	(0020,4000)	USER	CONDITIONAL	CONDITIONAL	ユーザー入力なしで元の値。入力した場合は入力値。	ユーザー入力の有無で異なる	
Instance Number	(0020,0013)	GENERATED	ALWAYS	ALWAYS	元の値に 10000 を足した値		
Burned InAnnotation	(0028,0301)	FIXED	ALWAYS	ALWAYS	NO		
Lossy Image Compression	(0028,2110)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Lossy Image Compression Ratio	(0028,2112)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Lossy Image Compression Method	(0028,2114)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

表 A-3-8 General Reference Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Derivation Description	(0008,2111)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Source Image Sequence	(0008,2112)	FIXED	ALWAYS	ALWAYS			
Referenced SOP Class UID	(0008,1150)	SRC_INSTANCE	SRC_COPY	SRC_COPY			
Referenced SOP Instance UID	(0008,1155)	SRC_INSTANCE	SRC_COPY	SRC_COPY			

表 A-3-9 SC Image Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Date Of Secondary Capture	(0018,1012)	GENERATED	ALWAYS	ALWAYS	DICOM を生成した日付		
Time Of Secondary Capture	(0018,1014)	GENERATED	ALWAYS	ALWAYS	DICOM を生成した時刻		

表 A-3-10 Image Pixel Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
Samples Per Pixel	(0028,0002)	GENERATED	ALWAYS	ALWAYS			
Photometric Interpretation	(0028,0004)	GENERATED	ALWAYS	ALWAYS			
Planar Configuration	(0028,0006)	FIXED	ALWAYS	ALWAYS	0		
Rows	(0028,0010)	GENERATED	ALWAYS	ALWAYS			
Columns	(0028,0011)	GENERATED	ALWAYS	ALWAYS			
Bits Allocated	(0028,0100)	FIXED	ALWAYS	ALWAYS	8		
Bits Stored	(0028,0101)	FIXED	ALWAYS	ALWAYS	8		
High Bit	(0028,0102)	FIXED	ALWAYS	ALWAYS	7		
Pixel Representation	(0028,0103)	FIXED	ALWAYS	ALWAYS	0		
Pixel Data	(7FE0,0010)	GENERATED	ALWAYS	ALWAYS			

表 A-3-11 SOP Common Module

Attribute Name	Tag	Source	Presence of Attribute	Presence of Value	Value	Conditions	Comments
SOP Class UID	(0008,0016)	FIXED	ALWAYS	ALWAYS	1.2.840.10008.5.1.4.1.1.7		Secondary Capture Image Storage
SOP Instance UID	(0008,0018)	GENERATED	ALWAYS	ALWAYS	1.2.392.200082.A.304.B		
Specific Character Set	(0008,0005)	SRC_INSTANCE / USER	CONDITIONAL	CONDITIONAL		使用する文字コードによる	

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Copyright (C) 2000-2024 GOODMAN Co., Ltd. All Rights Reserved.