

IHE Integration Statement			
Vendor	Product Name	Version	Date
Goodman Co., Ltd.	Goodnet	Version5	2024-3-29
This product implements all transactions required in the IHE Technical Framework to support the IHE profiles, actors and options listed below.			
Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Catheterization WorkFlow	Image Manager/Image Archive	None	
	Image Display	None	
Consistent Time	Time Client	None	
Retrieve ECG for Display	Display	None	
Echocardiography WorkFlow	Image Manager/Image Archive	None	
	Image Display	None	
Cardiology Evidence Documents	Image Manager/Image Archive	None	
	Image Display	None	
Portable Data for Imaging	Portable Media Creation	None	
	Portable Media Importer	None	
Patient Information Reconciliation	Image Manager/Image Archive	None	
Stress Testing WorkFlow	Image Manager/Image Archive	None	
	Image Display	None	
Scheduled WorkFlow	Image Manager/Image Archive	None	
	Image Display	None	
<b>Internet address for vendor's IHE information:</b>			
<a href="https://www.goodmankk.com/goodnet/download/index.html">https://www.goodmankk.com/goodnet/download/index.html</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>	<a href="https://www.goodmankk.com/goodnet/download/index.html">https://www.goodmankk.com/goodnet/download/index.html</a>		
<b>Links to general information about IHE worldwide: <a href="http://ihe.net/IHE_Worldwide">http://ihe.net/IHE_Worldwide</a></b>			