

IHE Integration Statement			
Vendor	Product Name	Version	Date
Goodman Co., Ltd.	Goodnet	Version7	2024-3-29
This product implements all transactions required in the IHE Technical Framework to support the IHE profiles, actors and options listed below.			
Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Catheterization WorkFlow	Image Manager/Image Archive	Intravascular Imaging	
	Image Display	Intravascular Imaging	
Consistent Time	Time Client	None	
Echocardiography WorkFlow	Image Manager/Image Archive	None	
	Image Display	None	
Patient Administration Management	Patient Demographics Consumer	None	
	Patient Encounter Consumer	None	
Portable Data for Imaging	Portable Media Creation	None	
	Portable Media Importer	None	
Patient Information Reconciliation	Image Manager/Image Archive	None	
Radiation Exposure Monitoring	Image Manager/Image Archive	None	
Scheduled WorkFlow.b	Image Manager/Image Archive	None	
	Image Display	None	
Endoscopy Image Archiving	Image Manager/Image Archive	None	
Internet address for vendor's IHE information: https://www.goodmankk.com/goodnet/download/index.html			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	https://www.goodmankk.com/goodnet/download/index.html		
Links to general information about IHE worldwide: http://ihe.net/IHE_Worldwide			